



Lenses and optical instruments

Laboratory course for students of the TUHH

Wladimir Banda-Barragán

Experimental procedure


Lens holder




Lenses




Lamp



Screen




Objects



Experimental procedure

1. Determining the focal length of a lens:



The focal lengths of two different converging lenses are to be determined by measuring the image distance **b** and object distance **g** for each lens.

Experimental procedure

1. Determining the focal length of a lens:





First a parallel beam of collimated light is produced using a lamp with the double condenser lens.

Next the object is positioned close to the condenser lens and a converging lens is used to produce a sharp image on the screen.

Experimental procedure


1. Determining the focal length of a lens:



The image distance and the object distance are measured from the middle of the thin lens.

Experimental procedure


2. Bessel method:



There are two different possible positions I and II in which a thin lens can produce a sharp image of an object at a distance d from the screen.

Experimental procedure

2. Bessel method:



For the lens in position I a sharp magnified image BI is produced on the screen.

For the lens in position II the object distance is equal to the image distance of the lens in position I and a sharp demagnified image BII is produced on the screen.

Experimental procedure

2. Bessel method:

Hence, at the two possible positions the image and object distances are both exchanged (i.e., $g_1 = b_2$ and $g_2 = b_1$).

$$g_1 + b_1 = d \quad (3a)$$



$$b_1 - g_1 = e \quad (3b)$$

➡ $b_1 = \frac{1}{2}(d + e) \quad \text{and} \quad g_1 = \frac{1}{2}(d - e)$

➡ $f = \frac{d^2 - e^2}{4d} \quad (4)$

Experimental procedure

2. Bessel method:




Thus, the focal length of a converging lens f_s can be determined from the values of d and e .

Experimental procedure

2. Bessel method:

This method can also be used to determine the focal length of a diverging lens f_z by measuring the focal length of a lens system f_{ges} consisting of a diverging lens and a converging lens:

$$\frac{1}{f_z} = \frac{1}{f_{ges}} - \frac{1}{f_s} \quad \text{oder} \quad f_z = \frac{f_{ges} \cdot f_s}{f_s - f_{ges}} \quad (5)$$



The power of the converging lens $1/f_s$ must be greater than that of the diverging lens $1/f_z$ to produce a real image.